

PPM1L antibody



Catalog Number: 114137

Product name

PPM1L antibody

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PPM1L Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue. Positive IHC detected in human kidney tissue, human heart tissue. Observed molecular weight by Western-blot: 41-45kd

Preparation

This antibody was obtained by immunization of PPM1L recombinant protein (Accession Number: NM_001317911). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:114137(PPM1L antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog No:114137(PPM1L Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human kidney tissue slide using Catalog

PPM1L antibody



Catalog Number: 114137

No:114137(PPM1L Antibody) at dilution of 1:50 (under 40x lens)