

# PPM1M antibody



Catalog Number: 114138

## Product name

PPM1M antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PPM1M Rabbit Polyclonal antibody. Positive WB detected in human lung tissue. Positive IP detected in mouse lung tissue. Observed molecular weight by Western-blot: 30 kDa

## Preparation

This antibody was obtained by immunization of PPM1M recombinant protein (Accession Number: NM\_001122870). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114138(PPM1M antibody) at dilution of 1:500



IP Result of anti-PPM1M (IP:Catalog No:114138, 4ug; Detection:Catalog No:114138 1:300) with mouse lung tissue lysate 4000ug.