

# PPM2C antibody

Catalog Number: 114139

## Product name

PPM2C antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PPM2C Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 53 kDa

## Preparation

This antibody was obtained by immunization of PPM2C recombinant protein (Accession Number: XM\_047421906). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:114139(PPM2C Antibody) at dilution of 1:300



IP Result of anti-PPM2C (IP:Catalog No:114139, 4ug; Detection:Catalog No:114139 1:300) with mouse brain tissue lysate 5200ug.