

PPOX antibody

Catalog Number: 114141

Product name

PPOX antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Immunogen

[Human PPOX Recombinant protein \(His tag\)](#)

Specificity

human, mouse, rat; other species not tested.

Antibody description

PPOX Rabbit Polyclonal antibody. Positive IHC detected in human kidney. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 51kd

Preparation

This antibody was obtained by immunization of PPOX recombinant protein (Accession Number: NM_000309). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:114141(PPOX antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human kidney using Catalog No:114141(PPOX antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:114141(PPOX

PPOX antibody



Catalog Number: 114141

antibody) at dilution of 1:50 (under 40x lens)