

PPWD1 antibody



Catalog Number: 114160

Product name

PPWD1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PPWD1 Rabbit Polyclonal antibody. Positive WB detected in Sp2/0 cells, HeLa cells, mouse brain tissue, mouse kidney tissue, mouse liver tissue, mouse testis tissue. Observed molecular weight by Western-blot: 65-66kd

Preparation

This antibody was obtained by immunization of PPWD1 recombinant protein (Accession Number: NM_015342). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Sp2/0 cells were subjected to SDS PAGE followed by western blot with Catalog No:114160(PPWD1 antibody) at dilution of 1:1000