

# PRCP antibody



Catalog Number: 114171

## Product name

PRCP antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PRCP Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, human liver tissue. Positive IP detected in L02 cells. Positive IHC detected in human lung cancer tissue. Observed molecular weight by Western-blot: 56 kDa

## Preparation

This antibody was obtained by immunization of PRCP recombinant protein (Accession Number: NM\_005040). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114171 (PRCP antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:114171 (PRCP antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded

# PRCP antibody



Catalog Number: 114171

human lung cancer using Catalog  
No:114171(PRCP antibody) at dilution of 1:50

IP Result of anti-PRCP (IP:Catalog No:114171,  
4ug; Detection:Catalog No:114171 1:300) with  
L02 cells lysate 1800ug.

