

# PRDM4 antibody



Catalog Number: 114174

---

## Product name

PRDM4 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

PRDM4 Rabbit Polyclonal antibody. Positive WB detected in mouse pancreas tissue, COLO 320 cells. Observed molecular weight by Western-blot: 88kd

## Preparation

This antibody was obtained by immunization of PRDM4 recombinant protein (Accession Number: NM\_012406). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:114174(PRDM4 antibody) at dilution of 1:500