

pregnancy zone protein antibody



Catalog Number: 114181

Product name

pregnancy zone protein antibody

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

pregnancy zone protein Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue, human brain tissue. Positive WB detected in human placenta tissue, human brain tissue. Observed molecular weight by Western-blot: 164kd

Preparation

This antibody was obtained by immunization of pregnancy zone protein recombinant protein (Accession Number: XM_024449111). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:114181(PZP antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human placenta using Catalog No:114181(PZP antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:114181(PZP antibody) at dilution of 1:50 (under 40x lens)