

# PRELID1 antibody



Catalog Number: 114183

## Product name

PRELID1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PRELID1 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human breast cancer tissue, human colon cancer tissue, human liver cancer tissue.

## Preparation

This antibody was obtained by immunization of PRELID1 recombinant protein (Accession Number: NM\_013237). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114183(PRELID1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:114183(PRELID1 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog

# PRELID1 antibody



Catalog Number: 114183

---

No:114183(PRELID1 Antibody) at dilution of 1:200 (under 40x lens).