

PREP antibody



Catalog Number: 114184

Product name

PREP antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PREP Rabbit Polyclonal antibody. Positive IHC detected in human gliomas. Positive WB detected in human lung tissue, Jurkat cells, mouse lung tissue. Observed molecular weight by Western-blot: 81kd

Preparation

This antibody was obtained by immunization of PREP recombinant protein (Accession Number: NM_002726). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



human lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:114184(PREP antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114184(PREP antibody) at dilution of 1:50 (under 10x lens)