

PRICKLE1 antibody

Catalog Number: 114190

Product name

PRICKLE1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PRICKLE1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue, Neuro-2a cells. Observed molecular weight by Western-blot: 94 kDa, 55-60 kDa

Preparation

This antibody was obtained by immunization of PRICKLE1 recombinant protein (Accession Number: XM_047428334). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114190(PRICKLE1 antibody) at dilution of



1:1000

IP Result of anti-PRICKLE1 (IP:Catalog No:114190, 3ug; Detection:Catalog No:114190 1:1000) with mouse brain tissue lysate 7000ug.