

# PRKAR1B antibody



Catalog Number: 114196

## Product name

PRKAR1B antibody

## Immunogen

[Human PRKAR1B Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PRKAR1B Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 45-50 kDa

## Preparation

This antibody was obtained by immunization of PRKAR1B recombinant protein (Accession Number: NM\_001164759). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB


## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114196(PRKAR1B Antibody) at dilution of



1:1000 mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114196(PRKAR1B Antibody) at dilution of 1:1000