

# PRMT3 antibody



Catalog Number: 114206

## Product name

PRMT3 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PRMT3 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 68 kDa

## Preparation

This antibody was obtained by immunization of PRMT3 recombinant protein (Accession Number: NM\_005788). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:114206 (PRMT3 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human brain using Catalog No:114206 (PRMT3 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human brain using Catalog No:114206 (PRMT3 antibody) at dilution of 1:100 (under 40x lens)