

# PRPF19 antibody



Catalog Number: 114227

## Product name

PRPF19 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PRPF19 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, human testis tissue, mouse kidney tissue, mouse liver tissue, PC-3 cells, rat kidney tissue. Observed molecular weight by Western-blot: 55kd

## Preparation

This antibody was obtained by immunization of PRPF19 recombinant protein (Accession Number: NM\_014502). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114227(PRPF19 antibody) at dilution of 1:500