

PRPF4 antibody



Catalog Number: 114231

Product name

PRPF4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PRPF4 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 58kd

Preparation

This antibody was obtained by immunization of PRPF4 recombinant protein (Accession Number: NM_004697). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:10-1:100

Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:114231 (PRPF4 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:114231 (PRPF4 antibody) at dilution of 1:10 (under 10x lens)