

PRPF40A antibody

Catalog Number: 114232

Product name

PRPF40A antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PRPF40A Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, mouse brain tissue. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 130-140 kDa

Preparation

This antibody was obtained by immunization of PRPF40A recombinant protein (Accession Number: BC027178). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114232 (PRPF40A Antibody) at dilution of



1:1000

IP Result of anti-PRPF40A (IP:Catalog No:114232, 4ug; Detection:Catalog No:114232 1:500) with HEK-293 cells lysate 3000ug.