

PSMF1 antibody



Catalog Number: 114250

Product name

PSMF1 antibody

Immunogen

[Human PSMF1 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PSMF1 Rabbit Polyclonal antibody. Positive WB detected in human kidney tissue, human heart tissue. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of PSMF1 recombinant protein (Accession Number: BC029836). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



human kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:114250(PSMF1 antibody) at dilution of 1:500