

PTBP2 antibody



Catalog Number: 114260

Product name

PTBP2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PTBP2 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, mouse brain tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 57-60 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_021190). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:114260(PTBP2 antibody) at dilution of 1:1000