

PTGR2 antibody



Catalog Number: 114273

Product name

Polyclonal

PTGR2 antibody

Ig Type

Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

Applications

Antibody description

ELISA, WB

PTGR2 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue. Observed molecular weight by Western-blot: 38.5kd

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of PTGR2 recombinant protein (Accession Number: NM_001146154). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:114273(PTGR2 antibody) at dilution of 1:500