

PRRG2 antibody



Catalog Number: 114277

Product name

Polyclonal

PRRG2 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

PRRG2 Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 16 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of PRRG2 recombinant protein (Accession Number: NM_001316335). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:114277(PRRG2 antibody) at dilution of 1:300