

# PSAP antibody



Catalog Number: 114291

## Product name

PSAP antibody

Recommended Dilution:

WB: 1:500-1:5000

## Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200


## Validations



## Antibody description

PSAP Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 60-70 kDa

PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114291(PSAP antibody) at dilution of 1:1500



## Preparation

This antibody was obtained by immunization of PSAP recombinant protein (Accession Number: NM\_002778). Purification method: Antigen affinity purified.

Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:114291(PSAP antibody) at dilution of 1:100 (under 10x lens)

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Immunohistochemical of paraffin-embedded human liver cancer using Catalog No:114291(PSAP antibody) at dilution of 1:100 (under 40x lens)

