

P-selectin antibody



Catalog Number: 114298

Product name

P-selectin antibody

Specificity

Human; other species not tested.

Antibody description

P-selectin Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue.

Preparation

This antibody was obtained by immunization of P-selectin recombinant protein (Accession Number: XM_005245438). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:114298(SELP antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:114298(SELP antibody) at dilution of 1:50 (under 40x lens)