

PTH-Specific antibody



Catalog Number: 114305

Product name

PTH-Specific antibody

Specificity

human; other species not tested.

Antibody description

PTH-Specific Rabbit Polyclonal antibody. Positive WB detected in Recombinant protein. Observed molecular weight by Western-blot: 16-19kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000315).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:114305(PTH-Specific antibody) at dilution of 1:500