

PTPMT1 antibody

Catalog Number: 114314

Product name

PTPMT1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PTPMT1 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer. Positive WB detected in mouse pancreas tissue. Observed molecular weight by Western-blot: 21-23kd

Preparation

This antibody was obtained by immunization of PTPMT1 recombinant protein (Accession Number: NM_175732). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse pancreas tissue were subjected to SDS PAGE followed by western blot with Catalog No:114314(PTPMT1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:114314(PTPMT1 antibody) at dilution of 1:50 (under 10x lens)