

PTPRJ antibody



Catalog Number: 114328

Product name

PTPRJ antibody

Specificity

Human; other species not tested.

Antibody description

PTPRJ Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, Jurkat cells, multi-cells/tissue. Observed molecular weight by Western-blot: 146-170kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002843). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications



ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



WB result of Catalog No:114328.

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:114328(PTPRJ antibody) at dilution of 1:500