

# PTPRR antibody



Catalog Number: 114329

---

## Product name

PTPRR antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

PTPRR Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, COLO 320 cells, HeLa cells. Observed molecular weight by Western-blot: 32kd, 47kd, 65kd

## Preparation

This antibody was obtained by immunization of PTPRR recombinant protein (Accession Number: NM\_130846). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114329(PTPRR antibody) at dilution of 1:500