

PURB antibody



Catalog Number: 114340

Product name

PURB antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PURB Rabbit Polyclonal antibody. Positive WB detected in C2C12 cell, mouse liver tissue. Observed molecular weight by Western-blot: 33-40 kDa

Preparation

This antibody was obtained by immunization of PURB recombinant protein (Accession Number: NM_033224). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



C2C12 cell were subjected to SDS PAGE followed by western blot with Catalog No:114340(PURB Antibody) at dilution of 1:600