

PXR antibody

Catalog Number: 114351

Product name

PXR antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PXR Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, A549 cells, MCF7 cells, mouse colon tissue. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 45-50kd

Preparation

This antibody was obtained by immunization of PXR recombinant protein (Accession Number: NM_003889). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:114351(NR1I2 antibody) at dilution of 1:200



IP Result of anti-PXR (IP:Catalog No:114351, 4ug; Detection:Catalog No:114351 1:500) with HepG2 cells lysate 3600ug.