

PYCR2 antibody



Catalog Number: 114356

Product name

PYCR2 antibody

Specificity

Human; other species not tested.

Antibody description

PYCR2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, L02 cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 27-31kd,20kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_013328). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114356(PYCR2 antibody) at dilution of 1:500



Immunofluorescent analysis of HeLa cells using Catalog No:114356(PYCR2 Antibody) at dilution of 1:50 and Rhodamine-Goat anti-Rabbit IgG