

PYGB antibody



Catalog Number: 114358

Product name

PYGB antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PYGB Rabbit Polyclonal antibody. Positive WB detected in A549 cells, HeLa cells, MCF7 cells. Positive IF detected in A549 cells. Observed molecular weight by Western-blot: 97-105kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002862).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:114358(PYGB antibody) at dilution of 1:500



Immunofluorescent analysis of A549 cells using Catalog No:114358(PYGB Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)