

PYK2 antibody

Catalog Number: 114363

Product name

PYK2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

PYK2 Rabbit Polyclonal antibody. Positive IHC detected in human medulloblastoma tissue, human kidney tissue. Positive WB detected in mouse brain tissue, Jurkat cells. Observed molecular weight by Western-blot: 116 kDa

Preparation

This antibody was obtained by immunization of PYK2 recombinant protein (Accession Number: XM_047421538). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114363(PTK2B antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human medulloblastoma using Catalog No:114363(PTK2B antibody) at dilution of 1:50 (under 10x lens)