

PSMB1 antibody

Catalog Number: 114379

Product name

PSMB1 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PSMB1 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, human brain tissue. Positive IF detected in HeLa cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 25kd

Preparation

This antibody was obtained by immunization of PSMB1 recombinant protein (Accession Number: NM_002793). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114379(PSMB1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:114379(PSMB1 antibody) at dilution of 1:50 (under 10x lens)



Immunofluorescent analysis of HeLa cells, using

PSMB1 antibody



Catalog Number: 114379

PSMB1 antibody Catalog No:114379 at 1:25
dilution and Rhodamine-labeled goat anti-rabbit

IgG (red). Blue pseudocolor = DAPI (fluorescent
DNA dye).