

PSMD13 antibody



Catalog Number: 114396

Product name

PSMD13 antibody

Specificity

Human; other species not tested.

Antibody description

PSMD13 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, Jurkat cells. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of PSMD13 recombinant protein (Accession Number: NM_002817). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114396(PSMD13 Antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells, using PSMD13 antibody Catalog No:114396 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).