

# RAB11FIP3 antibody



Catalog Number: 114409

## Product name

RAB11FIP3 antibody

## Immunogen

[Human RAB11FIP3 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RAB11FIP3 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 82 kDa

## Preparation

This antibody was obtained by immunization of RAB11FIP3 recombinant protein (Accession Number: NM\_014700). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114409(RAB11FIP3 Antibody) at dilution of 1:600