

RAB11FIP5 antibody



Catalog Number: 114410

Product name

RAB11FIP5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAB11FIP5 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HEK-293 cells, human heart tissue, mouse lung tissue. Observed molecular weight by Western-blot: 66-72 kDa

Preparation

This antibody was obtained by immunization of RAB11FIP5 recombinant protein (Accession Number: NM_015470). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:114410(RAB11FIP5 antibody) at dilution of 1:500