

RAB17 antibody

Catalog Number: 114414

Product name

RAB17 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAB17 Rabbit Polyclonal antibody. Positive IF detected in WIF-B cells.

Preparation

This antibody was obtained by immunization of RAB17 recombinant protein (Accession Number: NM_022449). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF


Dilutions

Recommended Dilution:

IF: N/A

Validations

rab17 antibody



IF result of RAB17 antibody (Catalog No:111366) with WIF-B cells overexpressed myc-tagged RAB17 by A.C. Striz and Dr. P.L. Tuma.