

# RAB1A antibody



Catalog Number: 114416

## Product name

RAB1A antibody

## Immunogen

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RAB1A Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human brain tissue, NIH/3T3 cells, rat brain tissue. Positive IHC detected in human pancreas cancer tissue, human liver tissue. Observed molecular weight by Western-blot: 23 kDa

## Preparation

This antibody was obtained by immunization of RAB1A recombinant protein (Accession Number: NM\_004161). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:114416(RAB1A antibody) at dilution of 1:50



(under 10x lens)

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114416(RAB1A Antibody) at dilution of 1:1000