

RAB1B antibody



Catalog Number: 114417

Product name

RAB1B antibody

Immunogen

[Human RAB1B Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

RAB1B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A375 cells, A431 cells. Positive IF detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 22-25 kDa

Preparation

This antibody was obtained by immunization of RAB1B recombinant protein (Accession Number: NM_030981). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114417(RAB1B antibody) at dilution of 1:500



Immunofluorescent analysis of NIH/3T3 cells using Catalog No:114417(RAB1B Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)