

# RAB1B antibody



Catalog Number: 114417

## Product name

RAB1B antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

RAB1B Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A375 cells, A431 cells. Positive IF detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 22-25 kDa

## Preparation

This antibody was obtained by immunization of RAB1B recombinant protein (Accession Number: NM\_030981). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114417(RAB1B antibody) at dilution of 1:500



Immunofluorescent analysis of NIH/3T3 cells using Catalog No:114417(RAB1B Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)