

RAB31 antibody



Catalog Number: 114429

Product name

RAB31 antibody

Immunogen

Specificity

human, mouse, rat; other species not tested.

Antibody description

RAB31 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 22kd

Preparation

This antibody was obtained by immunization of RAB31 recombinant protein (Accession Number: NM_006868). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:114429(RAB31 antibody) at dilution of 1:300