

# RAB31 antibody



Catalog Number: 114429

## Product name

Polyclonal

RAB31 antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

RAB31 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 22kd

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of RAB31 recombinant protein (Accession Number: NM\_006868). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:114429(RAB31 antibody) at dilution of 1:300