

# RAB40B antibody

Catalog Number: 114449

## Product name

RAB40B antibody

## Specificity

Human; other species not tested.

## Antibody description

RAB40B Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IHC detected in human gliomas tissue.

## Preparation

This antibody was obtained by immunization of RAB40B recombinant protein (Accession Number: NM\_006822). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114449(RAB40B antibody) at dilution of 1:50 (under 10x lens)



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114449(RAB40B Antibody) at dilution of 1:1000