

RAB41 antibody



Catalog Number: 114450

Product name

RAB41 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

RAB41 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, human brain tissue, human kidney tissue, mouse small intestine tissue, mouse uterus tissue. Observed molecular weight by Western-blot: 25kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001032726). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114450(RAB41 antibody) at dilution of 1:1000