

RAB4A antibody



Catalog Number: 114452

Product name

RAB4A antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAB4A Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue, HeLa cells, mouse brain tissue, rat brain tissue. Positive IF detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 26 kDa

Preparation

This antibody was obtained by immunization of RAB4A recombinant protein (Accession Number: NM_004578). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:20-1:200

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114452(RAB4A Antibody) at dilution of 1:300



Immunofluorescent analysis of (-20oc Ethanol) fixed NIH/3T3 cells using Catalog No:114452(RAB4A Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)