

# RAB5B antibody



Catalog Number: 114455

---

## Product name

RAB5B antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RAB5B Rabbit Polyclonal antibody. Positive WB detected in A431 cells, U-937 cells. Observed molecular weight by Western-blot: 25 kDa

## Preparation

This antibody was obtained by immunization of RAB5B recombinant protein (Accession Number: XM\_005269051). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:114455(RAB5B antibody) at dilution of 1:300