

RAB8B antibody



Catalog Number: 114462

Product name

RAB8B antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

RAB8B Rabbit Polyclonal antibody. Positive WB detected in mouse skin tissue, human brain tissue, mouse colon tissue. Observed molecular weight by Western-blot: 28kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_016530). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse skin tissue were subjected to SDS PAGE followed by western blot with Catalog No:114462(RAB8B antibody) at dilution of 1:1000