

RAB9A-Specific antibody



Catalog Number: 114464

Product name

RAB9A-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RAB9A-Specific Rabbit Polyclonal antibody.
Positive WB detected in HEK-293 cells, HeLa cells.
Positive IHC detected in human gliomas tissue.
Observed molecular weight by Western-blot: 23 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: XM_047442644).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114464(RAB9A-Specific antibody) at dilution of



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:114464(RAB9A-Specific Antibody) at dilution of 1:50