

# RARG antibody

Catalog Number: 114469

## Product name

RARG antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RARG Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in HeLa cells, MCF-7 cells. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of RARG recombinant protein (Accession Number: BC019098). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114469(RARG antibody) at dilution of 1:400



IP Result of anti-RARG (IP:Catalog No:114469, 4ug; Detection:Catalog No:114469 1:300) with MCF-7 cells lysate 1120ug.