

RASA3 antibody

Catalog Number: 114471

Product name

RASA3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RASA3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 100 kDa

Preparation

This antibody was obtained by immunization of RASA3 recombinant protein (Accession Number: NM_007368). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114471(RASA3 antibody) at dilution of 1:500