

RASGRP2 antibody



Catalog Number: 114478

Product name

RASGRP2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

RASGRP2 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, mouse liver tissue, mouse lung tissue. Observed molecular weight by Western-blot: 48kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_153819). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:114478(RASGRP2 antibody) at dilution of 1:300